



Ultrasound Telemedicine Made Simple

Share live scanning with anyone by SMS or email


One Button Telemedicine

All you need is the Clarius App to start a telemedicine session.



No 3rd party account or app needed

Simply send an SMS message or email from within the Clarius App to connect.



1. Start the Clarius App and your scanner

Once you're up and running, tap the Clarius Live icon.

2. Select a phone number or email

Use your device's address book or enter it manually.

3. Have recipient open the link on their device

On the other end, your recipient will receive a secure link that will instantly start the session.

View in a browser, on any device

Your recipients don't need the app to view sessions. Clarius Live sends a link that anyone can open in a web browser, on whatever device they're using.



First reactions

Check out what people had to say when we demoed Clarius Live at AIUM 2019 in Orlando, Florida.



Remote feedback

Get precise feedback. Recipients can point to areas of interest with their cursor, and the pointer will be highlighted on the scanning screen.

